


# Pulverizing Mill Bowl for LM2 and LM2000 Mills

This pulverizing mill bowl is designed for use with LM2 and LM2000 mills. It is suitable for processing sample masses ranging from 300 to 1600 grams and is available in standard steel and chrome options.



## ADDITIONAL IMAGES



## Product Overview

### High-Performance Pulverising Mill Bowl

This pulverising mill bowl features a professional-grade 'bowl and disc' style grinding head, making it a preferred choice for leading commercial assay laboratories. Engineered for durability, the unit is constructed from through-hardened steel to ensure an extended operational life under demanding industrial conditions. It is specifically designed for compatibility with LM2 and LM2000 mills, offering consistent particle size reduction and efficient material processing.

## Technical Specifications

Nominal Capacity	2000 cc
Recommended Sample Mass	300 to 1600 g
Total Bowl Weight	27.9 kg
Disc Weight	11.6 kg

## Dimensions

### Bowl Dimensions

**258 mm**  
Diameter

**93 mm**  
Height

## Compatibility & Materials

### Available Steel Types

Standard Steel • Chrome Steel

Compatible Mills	LM2, LM2000
------------------	-------------